



NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORIZATION FROM VERMEIREN GROUP

<b>INOVYS 2 EVO BACK FRAME</b>			
Scale:	1:5	Drafted	Verified
Date	20/04/2022	24/10/2025	24/10/2025
Initials	MK	PM	AS

---
Plan No.
E00035138/E1.1


Material			
No.	Plan	Article	Quantity
1	M8 Nut DIN985 Black	1996213	2
2	X001052	G066654	2
3	Washer M8 DIN126 Black	1996214	2
4	X003646	1999963	2
5	X004833	1745099	4
6	X004832	1744999	4
7	X004822	1744799	2
8	X004823	1744899	2
9	X004861	1906454	2
10	Y004353	F000411	2
11	Y017806	1756254 (40)	1
	Y017651	1756054 (45)	1
	Y017807	1756354 (50)	1
	YP06467	2734754 (35)	1
12	Y018198	1765462 (40)	1
	Y018199	1765562 (45)	1
	Y018200	1765662 (50)	1
	YP06468	2734862 (35)	1
13	Star Knob M8x20	1001580	2
14	Washer M8 PA6 Black	1918718	2
15	3,9x16 DIN7982	3994497	4
16	M6x25 ISO7380 Black	2994062	2
17	M6x25 DIN912 Black	2996887	4
18	Star Knob M6x20	1009932	2
19	Cap Ø22	1990648	4
20	Pin D8x40	1992510	1
21	X005513	1992600	2
22	---	1996806	2
23	M6x16 DIN7991 Black	2994064	4
24	M8x20 ISO7380 Black	1915944	2
25	L34	2995481 (35)	1
	L34	2995482 (40)	1
	L34	2995483 (45)	1
	L34	2995484 (50)	1

(40 - 45 - 50 - 35 = Seat width)

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

<h2 style="text-align: center;">INOVYS 2 EVO BACK FRAME</h2>				---
				Plan No. E00035138/E1.1
Scale: 1:5	Drafted	Verified	Confirmed	
Date	20/04/2022	24/10/2025	24/10/2025	
Initials	MK	PM	AS	