




Bill of Material			
No.	Plan	Article	Quantity
1	Side frame	E00050067	1
2	Cross	E00050163	1
3	Rear wheels	E00051594	1
4	Front wheels	E00051691	1
5	Backrest	E00049624	1
6	Sports brake	E00050162	1
	Push/Pull brake	E00050534	1
7	B06	E00048657	1
	B06 IN	EC00219	1
	*B06 / 90°	E00087278	1
	BA6	EC00305	1
	BK6	E00050417	1
	B08	E000981	1
	BK7/BL7	EP00499	1
8	B82(fender)	E00050783	1
	Bi5/Bi5A	E00050770	1
9	Upholstery	E00048650	1
	Upholstery options	A00100761	1
	L40	E00051931	1
10	B78	E00045457	1
11**	B90	E000734	1
12	B31	E00094427	1
13**	Tubular armrest	E00121056	1
14**	B75	E00129127	1

* - use only with front wheel`s size 4" or 5"

** - not depicted

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS
 PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN
 WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE,
 EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION
 FROM VERMEIREN GROUP.

V500 ACTIVE EXPLODED

Scale: 1:8	Drafted	Verified	Confirmed
Date	03/02/2023	16/06/2026	16/06/2026
Initials	AS	PM	AS

--
Plan No.
E00051729/F1.1
