


\* - used only with SE19

EP00489 / G


SE53 WIDE ARMREST / TURIOS				Plan No. EP00489 / G
Scale: 1:5	Drafted	Verified	Confirmed	
Date	12/02/2020	2/01/2024	2/01/2024	
Initials	DR	PM	MK	

Material				Material			
No.	Plan	Article	Quantity	No.	Plan	Article	Quantity
1	-	2606454	1	28*	M5X6	1997822	2
2	Y003626	1476554	1	29	Linx DLX REM 050	1909020	1
3	Y009175	1641654	3		Linx DLX REM 060	1910700	1
4	Y009176	1641754	1		Linx SE19 SE39	1912296	1
5	YP08065	2605954	2		Linx SE39	1912301	1
6A	YP08058	2605454 (L)	1		Linx SE19 ACT5	1912305	1
	YP08059	2605554 (R)	1		Linx LCD +JOYST	1913523	1
6B*	YP08060	2605654 (L)	1	30	M5x16	1907493	2
	YP08061	2605754 (R)	1	31*	Light controller	1915936	1
7*	YP08063	2605854	1	32	Spring plunger	1991432	1
8	YP09265	2642654	1	33	-	-	-
9	-	-	-	34	Armrest wide Soft Gray	1999993 (R)	1
10	X005565	1912287	1		Armrest wide Soft Gray	1999992 (L)	1
11	Armrest right	1913465	1		Armrest wide Soft Black	2998254 (R)	1
	Armrest left	1913464	1		Armrest wide Soft Black	2998253 (L)	1
12	XP00227	2990147	2	35*	M3x8	2990309	2
13	-	-	-	36	M6x45	2994058	2
14	Washer M6	1000973	8	37	Tube Cap 16x16	2994195	2
15	M5x12	1001489	2	38	XP03080	2996656	1
16	Tube Cap 20x20	1001612	2				
17	Nut M10	1001672	1				
18	Nut M5	1001679	6				
19	Nut M6	1001680	6				
20	M6x40	1003111	2				
21	Star Knob	1003658	1				
22	M5x20	1003907	4				
23	M6x40	1006111	4				
24	-	-	-				
25	M6x16	1008459	2				
26	Nut M6	1008492	2				
27	-	-	-				

Bill of Material 1913465/1913464 for EP00489			
No.	Plan	Article	Quantity
11a	X005567	1912286 (R)	1
	X005566	1912285 (L)	1
11b	X005517	2014999	1
11c	X005561	2021699	1
11d	M6x16	1008459	1
11e	Nut M6	1008492	1
11f	Rivet Ø4	1992805	2
11g	M5x16	3991136	2
11h	Washer M5	1906561	2
11i	Sticker	T001002	1

\* - used only with SE19

EP00489 / G

SE53 WIDE ARMREST / TURIOS				Plan No.
				EP00489 / G
Scale: 1:5	Drafted	Verified	Confirmed	
Date	12/02/2020	2/01/2024	2/01/2024	
Initials	DR	PM	MK	