





NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS
 PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN
 WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE,
 EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION
 FROM VERMEIREN GROUP.

Explosion Sagitta			
Scale:	1:7	Drafted	Verified
Date	31.07.2008	3.07.2026	3.07.2026
Initials	RM	PM	AS

Plan No.
E000481/O


Isometric view
Scale: 1:2



Front view
Scale: 1:5



Wheel / Kolko Ø 4"

Wheel / Kolko Ø 3"
SAGITTA LOW
SAGITTA LOW LONG

Wheel / Kolko Ø 4"
SAGITTA LOW
SAGITTA LONG

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS
PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN
WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE,
EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION
FROM VERMEIREN GROUP.


Explosion Sagitta

Scale: 1:7	Drafted	Verified	Confirmed
Date	31.07.2008	3.07.2026	3.07.2026
Initials	RM	PM	AS

Plan No.

E000481/O






NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS
PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN
WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE,
EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION
FROM VERMEIREN GROUP.

Explosion Sagitta			
Scale:	1:7	Drafted	Verified
Date	31.07.2008	3.07.2026	3.07.2026
Initials	RM	PM	AS

Plan No.

E000481/O





NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS
PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN
WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE,
EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION
FROM VERMEIREN GROUP.

Explosion Sagitta


Scale: 1:7	Drafted	Verified	Confirmed
Date	31.07.2008	3.07.2026	3.07.2026
Initials	RM	PM	AS

4/7

Plan No.

E000481/O






NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS
PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN
WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE,
EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION
FROM VERMEIREN GROUP.


Explosion Sagitta			
Scale:	1:7	Drafted	Verified
Date	31.07.2008	3.07.2026	3.07.2026
Initials	RM	PM	AS

Plan No.
E000481/O







0°




2°



4°



6°



8°

AXLES OF WEELS / OSIE KOL	
Lenght / Dlugosc L (mm)	Width / Szerokosc wozka (cm)
10	36
30	38
50	40
70	42
90	44
100	46

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS
PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN
WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE,
EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION
FROM VERMEIREN GROUP.

Explosion Sagitta			
Scale: 1:7	Drafted	Verified	Confirmed
Date	31.07.2008	3.07.2026	3.07.2026
Initials	RM	PM	AS

Plan No.
E000481/O


Material				Material				Material			
No./Nr	Plan/Nr rys.	Article/Artykul	Number/Ilosc	No./Nr.	Plan/Nr rys.	Article/Artykul	Number/Ilosc	No/Nr.	Plan/Nr rys.	Article/Artykul	Number/Ilosc
1	Frame	E000401	-	31	XP00392	GP00055	2	59	L	1992929	1
2	Footrest			32	D6	2990462	2	60	R	1992930	1
3	Backrest			33		1001557	2	61	YP01689	18766	1
4	YP01268	18534	1	34		1001554	2	*	YP05038*	28981*	1*
5	XP00353	1862554	1	35	M8	1001682	2	62	YP02932	28181	2
6	XP00352	1862654	2	36	M5x5	1991323	6	*	YP05037*	28980*	2*
7	XP00351	1862754	1	37		1900072	2	63	XP00366	2990331	4
8	---	---	---	38	YP01730	18772	1	64	YP01715	18767	1
	XP00454	2990349	2	39	YP01731	18897	1	*	YP05039*	28982*	1*
9	D4x8	2990307	2	40	XP00365	18553	4	65	XP00412	2990456	4
10		2990317	2	41	XP00357	2990300	4	66	D4X14	1001584	7
11	MBL628278	2990330(24")	2	42	Y000462	F000050	2	67	XP00422	2990457	4
	MBL618353	1991876(22")	2	43	X000273	G060600	2	68	M6x40	2990542	4
12	XP00385	18707(LOW)	2	44	YP01257	18533..	2	69	YP01328	18653	2
	XP00336	18628	2		YP01438	18631..	2	70		2990351	2
13	X000188	G0785	4		YP01439	18632..	2	71	YP01435	18624/18652	1
14		1003276	4		YP01440	18633..	2	72	M6x30	2990540	6
15	XP00338	18538	2	45		1001602	2	73	M6x25	1001431	2
16	XP00340	18548	2	46	GN717-6-AK	2990301	2	74	X001473	1002994	2
17	KNSTF 3"	2990568	2	47	M6x40 DIN 7991	2991312	2	75	4,8x16	1001534	2
	KNSTF 4"	2990348	2	48		2990164	12	76	---	1003299	0,9 M
	KNSTF 5"	2990328	2	49		1000973	2	** 77	BA6	E000410	---
18	XP00341	1854900	2	50	M6	1008492	8	** 78	PH05(Ø19)	EP00085	---
19	XP00339(4",5"-normal)	GP00044	4	51		1991857	2		PH05(Ø22)	EP00648	---
	XP00680 (4",5"-led;3",4"-trans)	2990917	4	52	D6xD18x1,6	1997328	2	** 79	B78	E000408	---
20		1001561	2	53	M6x25	2990538	2	** 80	PH03	EP00647	---
21	M12x1,25	1003369	2	54	M6x12	1005949	2				
22		1001564	4	55	XP00350	18550	2				
23	M6x16	2990537	4	56	MBL 674001	2990321	2				
24	XP00535	2990627	2	57	MBL 673074	2990319	1				
25	XP00345	2990310	2	58	MBL 673073	2990320	1				
26	XP00390	2990459	2								
27	M6x45 DIN 7991	2991311	2								
28	XP00405	GP00054	2								
29	XP00346	1853900	2								
30	XP00369	GP00045	4								

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

Explosion Sagitta

Scale: 1:7

Drafted

Verified

Confirmed

Date

31.07.2008

3.07.2026

3.07.2026

Initials

RM

PM

AS

Plan No.

E000481/O



* - option for France

** - not depicted