




No.	Plan	Article	Qty
1	E000634	1906314	1
2	P00083426	2691154	1
2**	P00122812	2923254	1
3	X004747	1906057	2
4	X004750	1906037	1
5	X004774	1906056	2
6	X004800	1740300	1
7	Y016850	1740132	1
8	M6 DIN985	1001680	1
9	---	1999935	1
10	M6 DIN985/ Black	2994063	2
11	M6x35	1006057	1
12	M6x25	1006998	2
13	D6,4 DIN125/Black	2994414	2
14	M6x35 DIN7984/Black	2994654	2

* - subparts of no.1 - E000634/1906314

** - Use with option amputee

B78 ANTI-TIPPING / V500 ACT/LIGHT

Scale: 1:5	Drafted	Verified	Confirmed
Date	2022-10-12	18/05/2026	18/05/2026
Initials	AZ	AS	AZ

Plan No.

E00045457/E1.1

