







Left view



Front view



Right view



1:20

2 Note: Serial number sticker must be fixed parallel to frame.

No.	Title	Art. nr.	Quantity
1	Sticker 35 logo Vermeiren	T001002	2
2	Sticker serial number	T001097	1
3	Sticker name vermeiren D2003990609		2

STICKER CONFIGURATION D200 30°

Scale: 2:25	Drafted	Verified	Confirmed
Date	20/05/2021	16/07/2026	16/07/2026
Initials	AZ	MP	MP

Plan Number

EC00504/B



EC00504/B