


use grease inside sleeve

Tighten nuts with torque wrench 5 Nm (W600135)

use with grease

Detail A

Without lights



NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

REAR LOWER FRAME TIMIX

Scale: 1:5	Drafted	Verified	Confirmed
Date	22.06.2012	23/06/2026	23/06/2026
Initials	JG	KA	AZ

Sheet 1/2

2991641

Plan No.

E000816/M



Material							
No.	Plan	Article	Quantity	No.	Plan	Article	Quantity
1	X001016	1001682	1	21	---	1907338	2
2	X001119	1001554	2	22	M12 ISO7379	1907488	1
4	X004951	1908263	1	23	---	1908027	1
5	X004967	1908262	2	24	---	1991679	2
6	X004976	1908264	2	25	M12x110 DIN912	1907346	2
7	X004999	1908130	1	26	Black Line Le.	1915166	1
8A	Y017903	17691..	1		Black R-fill Le.	1915167	1
8B*	Y019739	17864..	1	28	M6 DIN125A	1000973	2
9	Y019289	17833..	1	29*	---	1992835	4
10	Y019286	17832..	2	30	Y021635	2013454	2
11	---	1001409	2	31*	---	1001611	2
12	M4x20 DIN84	1001466	4	32	3,9X9,5 DIN7981	1001530	2
13	M12 DIN125B	1001561	5	33	X005589	1912549	2
14	M6 DIN125A	1001564	4	34	Black Line Ri.	1915346	1
15	M4 DIN985	1001678	4		Black R-fill Ri.	1915347	1
16	M6 DIN985	1001680	4				
17	M12 DIN985	1001684	5				
18	M6x45 ISO7380	1991017	2				
19	M6x50 DIN912	1007970	2				
20	M8x60 ISO7380	1900809	1				

* Only without lights

Sheet 2/2

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS
 PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN
 WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE,
 EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION
 FROM VERMEIREN GROUP.

2991641

REAR LOWER FRAME TIMIX

Plan No.

E000816/M

Scale: 1:5	Drafted	Verified	Confirmed
Date	22.06.2012	23/06/2026	23/06/2026
Initials	JG	KA	AZ

