



Material Inovy's Evo			
No.	Plan	Article	Quantity
1	Y011362	1658154 (40)	1
	Y011932	1663554 (45)	1
	Y011936	1663654 (50)	1
	YP06990	2743654 (35)	1
2	Y011922	1658354 (R)	1
3	Y011422	1658254 (L)	1
4	M8 DIN126	1996214	10
5	M8 DIN 985	1996213	6
6	M12 (PA6.6)	1916509	4
7	M8 (PA6.6)	2994415	4
8	M8x60	2995491	2
9	M8x55	1919452	2
10	Nut M8	1998994	2

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

Sheet 1/2

---

### INOVYS PIVOT FRAME

Plan No.

E000383/J

Scale: 1:1	Drafted	Verified	Confirmed
Date	5/07/2007	24/10/2025	24/10/2025
Initials	SD	PM	AS



Material Inovys II			
No.	Plan	Article	Quantity
1	Y011362	1658154 (40)	1
	Y011932	1663554 (45)	1
	Y011936	1663654 (50)	1
	YP06990	2743654 (35)	1
2	Y011922	1658354 (R)	1
3	Y011422	1658254 (L)	1
4	M8 DIN126	1001554	10
5	M8 DIN 985	1001682	6
6	M12 (PA6.6)	1916509	4
7	M8 (PA6.6)	1003199	4
8	M8x60	1900809	2
9	M8x55	1995854	2
10	Nut M8	1915942	2

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

Sheet 2/2

---

Plan No.

E000383/J

INOVYS PIVOT FRAME

Scale:	1:1	Drafted	Verified	Confirmed
Date	5/07/2007	24/10/2025	24/10/2025	
Initials	SD	PM	AS	

