


EXPLODED VIEW 275 BIS				Plan No.
				1/2 E000114/M
Scale: 1:8	Drafted	Verified	Confirmed	
Date	12/03/2004	29/11/2021	29/11/2021	
Init.	BS	AZ	TK	

Material			
No. / Nr.	Plan / Rysunek	Article / Artikel	Number / Aantal
1	YP09307DP	26435...	1
2	YP00372	18117...	1
3	Y004067	15072...	1
4	Y004072	15073...	2
5	YP09302	2643404	2
6	---	1611704	1
7	---	1900183	4
8	M12x60	1993419	2
9	---	1991447	2
10	M8x76 DIN912/8.8	1998993	2
11	M10x80 DIN931/8.8	1001457	2
12	M10 DIN985	1001683	2
13	M6x75 DIN912/8.8	1993413	2
14	X002306	1502304	2
15	X002443	1517104	2
16	Ø10x1x39.5	1517200	2
17	M6 DIN125 A	1001564	4
18	M6 DIN985	1001680	12
19	M6x40	1006111	2
20	M6x60 DIN912/8.8	1005957	4
21	PVC M6 DIN125 A	1000973	4
22	XP00099	2990067	4
23	M10 DIN125 A	1001558	6
24	M6x16	1008459	2
25	Ø4x30	1001545	4
26	X002305	1993120	2
27	Øx20	2990326	2
28	---	1909590	2
29	---	1908896	2
30	---	1908909	2
31	/	/	/

Material			
No. / Nr.	Plan / Rysunek	Article / Artykul	Number / Iloscl
32	/	/	/
33	/	/	/
34*	---	1908105	4
35	---	1909113	2
36	---	1993417	1
37	X005513	1992600	2
38	/	/	
39	---	1909576	2
40	L=960 mm 125 cm	1003279 1004094	2 2
41	XP00071	2990040	2
42	Øx12	1990816	4
43	---	1994587	1
44	X002335	1993411	1
45	Y004199	1511304	1
46	Y004200	1511404	1
47	X002170	1992792	2
48	XP00108	GP006...	2
49	M12 DIN127 A	1001572	2
50	M6x23	1003558	2
51	M5x10 DIN916/45H	1001485	2
52*	---	T001060	2
53*	---	T001061	2
54*	LOGO VERMEIREN	T001002	2
55	M8	1998994	2

\* - Not depicted

EXPLODED VIEW 275 BIS				Plan No.
				2/2 E000114/M
Scale: 1:8	Drafted	Verified	Confirmed	
Date	12/03/2004	29/11/2021	29/11/2021	
Init.	BS	AZ	TK	